

CYP2A6/7 rabbit pAb

Cat No.: ES4926

For research use only

Overview

Product Name CYP2A6/7 rabbit pAb

Host species Rabbit
Applications IHC;IF;ELISA

Species Cross-Reactivity Human; Rat; Mouse;

Recommended dilutions Immunohistochemistry: 1/100 - 1/300. ELISA:

1/40000. Not yet tested in other applications.

Immunogen Synthesized peptide derived from CYP2A6/7 . at AA

range: 310-390

Specificity CYP2A6/7 Polyclonal Antibody detects endogenous

levels of CYP2A6/7 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Cytochrome P450 2A7

Gene Name CYP2A7

Cellular localization Endoplasmic reticulum membrane; Peripheral

membrane protein. Microsome membrane;

Peripheral membrane protein.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

Observed band

Human Gene ID 1549 Human Swiss-Prot Number P20853

Alternative Names CYP2A7; Cytochrome P450 2A7; CYPIIA7;

Cytochrome P450 IIA4

Background This gene encodes a member of the cytochrome

P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This protein



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com



localizes to the endoplasmic reticulum; its substrate has not yet been determined. This gene, which produces two transcript variants, is part of a large cluster of cytochrome P450 genes from the CYP2A, CYP2B and CYP2F subfamilies on chromosome 19q. [provided by RefSeq, Jul 2008],



Immunohistochemical analysis of paraffin-embedded human Colon cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).



+86-27-59760950