



ELK Biotechnology

CYP2D6 Rabbit pAb
Catalog Number: ES4938
For research use only.

Overview

Product Name	CYP2D6 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/40000
Immunogen	Synthesized peptide derived from the Internal region of human CYP2D6. at AA range; 210-290
Specificity	CYP2D6 Polyclonal Antibody detects endogenous levels of CYP2D6 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	CYP2D6
Gene Name	CYP2D6
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	55kDa
Human Gene ID	1565
Human Swiss-Prot Number	P10635
Alternative Names	CYP2D6; CYP2DL1; Cytochrome P450 2D6; CYP1ID6; Cytochrome P450-DB1; Debrisoquine 4-hydroxylase