

Nkx-6.3 rabbit pAb

Cat No.:ES4945

For research use only

Overview

Product Name	Nkx-6.3 rabbit pAb
Host species	Rabbit
Applications	WB;IF;ELISA
Species Cross-Reactivity	Human; Mouse; Monkey
Recommended dilutions	Western Blot: 1/500 - 1/2000. Immunofluorescence:
	1/200 - 1/1000. ELISA: 1/10000. Not yet tested in
	other applications.
Immunogen	The antiserum was produced against synthesized
	peptide derived from human NKX6.3. AA
	range:161-210
Specificity	Nkx-6.3 Polyclonal Antibody detects endogenous
	levels of Nkx-6.3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Homeobox protein Nkx-6.3
Gene Name	NKX6-3
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	34kD
Human Gene ID	157848
Human Swiss-Prot Number	A6NJ46
Alternative Names	NKX6-3; Homeobox protein Nkx-6.3
Background	The NKX family of homeodomain proteins controls
	numerous developmental processes. Members of
	the NKX6 subfamily, including NKX6-3, are involved
	in development of the central nervous system (CNS),
	gastrointestinal tract, and pancreas (Alanentalo et
	al., 2006 [PubMed 16326147]).[supplied by OMIM,




+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C







Mar 2008],

Western Blot analysis of K562 cells using Nkx-6.3 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003,Inventbiotech,MN,USA).

Immunofluorescence analysis of HUVEC cells, using NKX6.3 Antibody. The picture on the right is blocked with the synthesized peptide.





Western blot analysis of lysates from K562 and COS7 cells, using NKX6.3 Antibody. The lane on the right is blocked with the synthesized peptide.



ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C