



ELK Biotechnology

CD55 Rabbit pAb
Catalog Number: ES4957
For research use only.

Overview

Product Name	CD55 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/20000
Immunogen	Synthesized peptide derived from the Internal region of human CD55. at AA range; 70-150
Specificity	CD55 Polyclonal Antibody detects endogenous levels of CD55 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	CD55
Gene Name	CD55
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	41kDa
Human Gene ID	1604
Human Swiss-Prot Number	P08174
Alternative Names	CD55; CR; DAF; Complement decay-accelerating factor; CD antigen CD55