

DDX3 rabbit pAb

Cat No.:ES4985

For research use only

Overview

Product Name	DDX3 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not
	yet tested in other applications.
Immunogen	The antiserum was produced against synthesized
	peptide derived from N-ternal human DDX3. AA
	range:14-63
Specificity	DDX3 Polyclonal Antibody detects endogenous
	levels of DDX3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	ATP-dependent RNA helicase DDX3X/Y
Gene Name	DDX3X/DDX3Y
Cellular localization	Cell membrane . Nucleus . Cytoplasm . Cytoplasm,
	Stress granule . Inflammasome . Cell projection,
	lamellipodium . Cytoplasm, cytoskeleton,
	microtubule organizing center, centrosome.
	Shuttles between the nucleus and the cytosol
	(PubMed:15507209, PubMed:18
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	73kD
Human Gene ID	1654/8653
Human Swiss-Prot Number	000571/015523
Alternative Names	DDX3X; DBX; DDX3; ATP-dependent RNA helicase
	DDX3X; DEAD box protein 3; X-chromosomal;
	DEAD box, X isoform; Helicase-like protein 2; HLP2;



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com


23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



Background

DDX3Y; DBY; ATP-dependent RNA helicase DDX3Y; DEAD box protein 3, Y-chromosomal The protein encoded by this gene is a member of the large DEAD-box protein family, that is defined by the presence of the conserved Asp-Glu-Ala-Asp (DEAD) motif, and has ATP-dependent RNA helicase activity. This protein has been reported to display a high level of RNA-independent ATPase activity, and unlike most DEAD-box helicases, the ATPase activity is thought to be stimulated by both RNA and DNA. This protein has multiple conserved domains and is thought to play roles in both the nucleus and cytoplasm. Nuclear roles include transcriptional regulation, mRNP assembly, pre-mRNA splicing, and mRNA export. In the cytoplasm, this protein is thought to be involved in translation, cellular signaling, and viral replication. Misregulation of this gene has been implicated in tumorigenesis. This gene has a paralog located in the nonrecombining region of the Y chromosome. Pseudogenes sharing similarit

Western Blot analysis of Jurkat cells using DDX3 Polyclonal Antibody diluted at 1:2000






www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



138 100 70 55 40 35 25	L929 MCF7 HELA	Western Blot analysis of L929 MCF7 HELA cells using Antibody diluted at 2000. Secondary antibody(catalog#:RS0002) was diluted at 1:20000
15		



Western blot analysis of lysate from Jurkat cells, using DDX3 antibody.



+86-27-59760950 ELKbio@ELKbiotech.com

com www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C