



## ELK Biotechnology

**CRMP-2 Rabbit pAb**  
**Catalog Number: ES5023**  
**For research use only.**

### Overview

---

Product Name	CRMP-2 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	<b>Western Blot:1/500-1/2000</b> <b>Immunohistochemistry:1/100-1/300</b> <b>ELISA:1/20000</b>
Immunogen	Synthesized peptide derived from human CRMP-2 around the non-phosphorylation site of T509.
Specificity	CRMP-2 Polyclonal Antibody detects endogenous levels of CRMP-2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	CRMP-2
Gene Name	DPYSL2
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	<b>50kDa</b>
Human Gene ID	1808
Human Swiss-Prot Number	Q16555
Alternative Names	DPYSL2; CRMP2; ULIP2; Dihydropyrimidinase-related protein 2; DRP-2; Collapsin response mediator protein 2; CRMP-2; N2A3; Unc-33-like phosphoprotein 2; ULIP-2