



ELK Biotechnology

Elk-1 Rabbit pAb
Catalog Number: ES5097
For research use only.

Overview

Product Name	Elk-1 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,IP,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 Immunoprecipitation:2-5ug/mg lysate ELISA:1/20000
Immunogen	Synthesized peptide derived from human Elk-1 around the non-phosphorylation site of S389.
Specificity	Elk-1 Polyclonal Antibody detects endogenous levels of Elk-1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Elk-1
Gene Name	ELK1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	45kDa
Human Gene ID	2002
Human Swiss-Prot Number	P19419
Alternative Names	ELK1; ETS domain-containing protein Elk-1