



# ELK Biotechnology

**Neu Rabbit pAb**  
**Catalog Number: ES5144**  
**For research use only.**

## Overview

---

<b>Product Name</b>	Neu polyclonal antibody
<b>Host species</b>	Rabbit
<b>Applications</b>	<b>WB,ELISA</b>
<b>Species Cross-Reactivity</b>	<b>Human,Mouse,Rat</b>
<b>Recommended dilutions</b>	<b>Western Blot:1/500-1/2000</b> <b>ELISA:1/10000</b>
<b>Immunogen</b>	Synthesized peptide derived from human Neu around the non-phosphorylation site of Y1112.
<b>Specificity</b>	Neu Polyclonal Antibody detects endogenous levels of Neu protein.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Protein Name</b>	Neu
<b>Gene Name</b>	ERBB2
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	<b>138kDa</b>
<b>Human Gene ID</b>	2064
<b>Human Swiss-Prot Number</b>	P04626
<b>Alternative Names</b>	ERBB2; HER2; MLN19; NEU; NGL; Receptor tyrosine-protein kinase erbB-2; Metastatic lymph node gene 19 protein; MLN 19; Proto-oncogene Neu; Proto-oncogene c-ErbB-2; Tyrosine kinase-type cell surface receptor HER2; p185erbB2; CD antigen CD340