



ELK Biotechnology

ERF Rabbit pAb
Catalog Number: ES5169
For research use only.

Overview

Product Name	ERF polyclonal antibody
Host species	Rabbit
Applications	IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse
Recommended dilutions	Immunohistochemistry:1/100-1/300 ELISA:1/40000
Immunogen	Synthesized peptide derived from the C-terminal region of human ERF. at AA range; 440-520
Specificity	ERF Polyclonal Antibody detects endogenous levels of ERF protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	ERF
Gene Name	ERF
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	2077
Human Swiss-Prot Number	P50548
Alternative Names	ERF; ETS domain-containing transcription factor ERF; Ets2 repressor factor; PE-2