



ELK Biotechnology

Fhit Rabbit pAb
Catalog Number: ES5254
For research use only.

Overview

Product Name	Fhit polyclonal antibody
Host species	Rabbit
Applications	IHC-p,WB,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/40000
Immunogen	Synthesized peptide derived from human Fhit around the non-phosphorylation site of Y114.
Specificity	Fhit Polyclonal Antibody detects endogenous levels of Fhit protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Fhit
Gene Name	FHIT
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	2272
Human Swiss-Prot Number	P49789
Alternative Names	FHIT; Bis(5'-adenosyl)-triphosphatase; AP3A hydrolase; AP3Aase; Diadenosine 5'; 5'''-P1,P3-triphosphate hydrolase; Dinucleosidetriphosphatase; Fragile histidine triad protein