



ELK Biotechnology

Foxl1 Rabbit pAb
Catalog Number: ES5283
For research use only.

Overview

Product Name	Foxl1 polyclonal antibody
Host species	Rabbit
Applications	IHC-p,WB,ELISA
Species Cross-Reactivity	Human,Mouse
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/5000
Immunogen	Synthesized peptide derived from the Internal region of human Foxl1. at AA range; 160-240
Specificity	Foxl1 Polyclonal Antibody detects endogenous levels of Foxl1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Foxl1
Gene Name	FOXI1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	2299
Human Swiss-Prot Number	Q12951
Alternative Names	FOXI1; FKHL10; FREAC6; Forkhead box protein I1; Forkhead-related protein FKHL10; Forkhead-related transcription factor 6; FREAC-6; Hepatocyte nuclear factor 3 forkhead homolog 3; HFH-3; HNF-3/fork-head homolog 3