



# ELK Biotechnology

**M-RIP Rabbit pAb**  
**Catalog Number: ES5313**  
**For research use only.**

## Overview

---

<b>Product Name</b>	M-RIP polyclonal antibody
<b>Host species</b>	Rabbit
<b>Applications</b>	<b>IHC-p,WB,ELISA</b>
<b>Species Cross-Reactivity</b>	<b>Human,Mouse,Rat</b>
<b>Recommended dilutions</b>	<b>Western Blot:1/500-1/2000</b> <b>Immunohistochemistry:1/100-1/300</b> <b>ELISA:1/40000</b>
<b>Immunogen</b>	Synthesized peptide derived from the Internal region of human M-RIP. at AA range; 700-780
<b>Specificity</b>	M-RIP Polyclonal Antibody detects endogenous levels of M-RIP protein.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Protein Name</b>	M-RIP
<b>Gene Name</b>	MPRIP
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	<b>110kDa</b>
<b>Human Gene ID</b>	23164
<b>Human Swiss-Prot Number</b>	Q6WCQ1
<b>Alternative Names</b>	MPRIP; KIAA0864; MRIP; RHOIP3; Myosin phosphatase Rho-interacting protein; M-RIP; Rho-interacting protein 3; RIP3; p116Rip