



ELK Biotechnology

HP1 α Rabbit pAb
Catalog Number: ES5347
For research use only.

Overview

Product Name	HP1 α polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/5000
Immunogen	Synthesized peptide derived from human HP1 α around the non-phosphorylation site of S92.
Specificity	HP1 α Polyclonal Antibody detects endogenous levels of HP1 α protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	HP1 α
Gene Name	CBX5
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	24kDa
Human Gene ID	23468
Human Swiss-Prot Number	P45973
Alternative Names	CBX5; HP1A; Chromobox protein homolog 5; Antigen p25; Heterochromatin protein 1 homolog alpha; HP1 alpha