



ELK Biotechnology

Encephalopsin Rabbit pAb
Catalog Number: ES5367
For research use only.

Overview

Product Name	Encephalopsin polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Cross-Reactivity	Human,Mouse
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 Immunofluorescence:1/200-1/1000 ELISA:1/10000
Immunogen	Synthesized peptide derived from the Internal region of human Encephalopsin. at AA range; 130-210
Specificity	Encephalopsin Polyclonal Antibody detects endogenous levels of Encephalopsin protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Encephalopsin
Gene Name	OPN3
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	48kDa
Human Gene ID	23596
Human Swiss-Prot Number	Q9H1Y3
Alternative Names	OPN3; ECPN; Opsin-3; Encephalopsin; Panopsin