



ELK Biotechnology

FSHR Rabbit pAb
Catalog Number: ES5399
For research use only.

Overview

Product Name	FSHR polyclonal antibody
Host species	Rabbit
Applications	WB,IF,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 Immunofluorescence:1/200-1/1000 ELISA:1/5000
Immunogen	Synthesized peptide derived from the Internal region of human FSHR. at AA range; 250-330
Specificity	FSHR Polyclonal Antibody detects endogenous levels of FSHR protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	FSHR
Gene Name	FSHR
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	78kDa
Human Gene ID	2492
Human Swiss-Prot Number	P23945
Alternative Names	FSHR; LGR1; Follicle-stimulating hormone receptor; FSH-R; Follitropin receptor