



ELK Biotechnology

Gab 1 Rabbit pAb
Catalog Number: ES5424
For research use only.

Overview

Product Name	Gab 1 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 Immunofluorescence:1/200-1/1000 ELISA:1/10000
Immunogen	Synthesized peptide derived from human Gab 1 around the non-phosphorylation site of Y659.
Specificity	Gab 1 Polyclonal Antibody detects endogenous levels of Gab 1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Gab 1
Gene Name	GAB1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	76kDa
Human Gene ID	2549
Human Swiss-Prot Number	Q13480
Alternative Names	GAB1; GRB2-associated-binding protein 1; GRB2-associated binder 1; Growth factor receptor bound protein 2-associated protein 1