



ELK Biotechnology

Pdcd-4 Rabbit pAb
Catalog Number: ES5551
For research use only.

Overview

Product Name	Pdcd-4 polyclonal antibody
Host species	Rabbit
Applications	IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Immunohistochemistry:1/100-1/300 ELISA:1/5000
Immunogen	Synthesized peptide derived from human Pdcd-4 around the non-phosphorylation site of S67.
Specificity	Pdcd-4 Polyclonal Antibody detects endogenous levels of Pdcd-4 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Pdcd-4
Gene Name	PDCD4
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	27250
Human Swiss-Prot Number	Q53EL6
Alternative Names	PDCD4; H731; Programmed cell death protein 4; Neoplastic transformation inhibitor protein; Nuclear antigen H731-like; Protein 197/15a