



ELK Biotechnology

Ang-1 Rabbit pAb
Catalog Number: ES5590
For research use only.

Overview

Product Name	Ang-1 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/10000
Immunogen	Synthesized peptide derived from the Internal region of human Ang-1.
Specificity	Ang-1 Polyclonal Antibody detects endogenous levels of Ang-1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Ang-1
Gene Name	ANGPT1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	60kDa
Human Gene ID	284
Human Swiss-Prot Number	Q15389
Alternative Names	ANGPT1; KIAA0003; Angiopoietin-1; ANG-1