



ELK Biotechnology

GPR17 Rabbit pAb
Catalog Number: ES5591
For research use only.

Overview

Product Name	GPR17 polyclonal antibody
Host species	Rabbit
Applications	WB,IF,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 Immunofluorescence:1/200-1/1000 ELISA:1/5000
Immunogen	Synthesized peptide derived from the Internal region of human GPR17. at AA range; 190-270
Specificity	GPR17 Polyclonal Antibody detects endogenous levels of GPR17 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	GPR17
Gene Name	GPR17
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	40kDa
Human Gene ID	2840
Human Swiss-Prot Number	Q13304
Alternative Names	GPR17; Uracil nucleotide/cysteinyl leukotriene receptor; UDP/CysLT receptor; G-protein coupled receptor 17; P2Y-like receptor; R12