



ELK Biotechnology

GPR150 Rabbit pAb
Catalog Number: ES5611
For research use only.

Overview

Product Name	GPR150 polyclonal antibody
Host species	Rabbit
Applications	IF,ELISA
Species Cross-Reactivity	Human,Mouse
Recommended dilutions	Immunofluorescence:1/200-1/1000 ELISA:1/5000
Immunogen	Synthesized peptide derived from the Internal region of human GPR150. at AA range; 30-110
Specificity	GPR150 Polyclonal Antibody detects endogenous levels of GPR150 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	GPR150
Gene Name	GPR150
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	285601
Human Swiss-Prot Number	Q8NGU9
Alternative Names	GPR150; Probable G-protein coupled receptor 150