

GPS2 rabbit pAb

Cat No.:ES5628

For research use only

Overview

Product Name	GPS2 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not
	yet tested in other applications.
Immunogen	The antiserum was produced against synthesized
<u> </u>	peptide derived from human GPS2. AA range:11-60
Specificity	GPS2 Polyclonal Antibody detects endogenous levels
	of GPS2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	G protein pathway suppressor 2
Gene Name	GPS2
Cellular localization	Nucleus . Mitochondrion . Cytoplasm, cytosol .
	Sumoylation regulates the subcellular location
	(PubMed:24943844). Relocates from the
	mitochondria to the nucleus following
	desumoylation, leading to mediate mitochondrial
	stress response (By similarity).
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	40kD
Human Gene ID	2874
Human Swiss-Prot Number	Q13227
Alternative Names	GPS2; G protein pathway suppressor 2; GPS-2
Background	This gene encodes a protein involved in G
	protein-mitogen-activated protein kinase (MAPK)
	signaling cascades. When overexpressed in



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



mammalian cells, this gene could potently suppress a RAS- and MAPK-mediated signal and interfere with JNK activity, suggesting that the function of this gene may be signal repression. The encoded protein is an integral subunit of the NCOR1-HDAC3 (nuclear receptor corepressor 1-histone deacetylase 3) complex, and it was shown that the complex inhibits JNK activation through this subunit and thus could potentially provide an alternative mechanism for hormone-mediated antagonism of AP1 (activator protein 1) function. [provided by RefSeq, Jul 2008],

Western Blot analysis of RAW264.7 cells using GPS2 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003,Inventbiotech,MN,USA).



RAW264.7

GPS2

138---

70----

55-

40---35---

> Western blot analysis of lysates from RAW264.7 cells, using GPS2 Antibody. The lane on the right is blocked with the synthesized peptide.



www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C