



ELK Biotechnology

GR Rabbit pAb
Catalog Number: ES5667
For research use only.

Overview

Product Name	GR polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/10000
Immunogen	Synthesized peptide derived from human GR around the non-phosphorylation site of S226.
Specificity	GR Polyclonal Antibody detects endogenous levels of GR protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	GR
Gene Name	NR3C1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	90kDa
Human Gene ID	2908
Human Swiss-Prot Number	P04150
Alternative Names	NR3C1; GRL; Glucocorticoid receptor; GR; Nuclear receptor subfamily 3 group C member 1