



ELK Biotechnology

Rac GAP1 Rabbit pAb
Catalog Number: ES5677
For research use only.

Overview

Product Name	Rac GAP1 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Cross-Reactivity	Human,Mouse
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 Immunofluorescence:1/200-1/1000 ELISA:1/20000
Immunogen	Synthesized peptide derived from the N-terminal region of human Rac GAP1. at AA range; 30-110
Specificity	Rac GAP1 Polyclonal Antibody detects endogenous levels of Rac GAP1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Rac GAP1
Gene Name	RACGAP1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	70kDa
Human Gene ID	29127
Human Swiss-Prot Number	Q9H0H5
Alternative Names	RACGAP1; KIAA1478; MGCRACGAP; Rac GTPase-activating protein 1; Male germ cell RacGap; MgcRacGAP; Protein CYK4 homolog; CYK4; HsCYK-4