



ELK Biotechnology

eEF2K Rabbit pAb
Catalog Number: ES5697
For research use only.

Overview

Product Name	eEF2K polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/40000
Immunogen	Synthesized peptide derived from human eEF2K around the non-phosphorylation site of S359.
Specificity	eEF2K Polyclonal Antibody detects endogenous levels of eEF2K protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	eEF2K
Gene Name	EEF2K
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	100kDa
Human Gene ID	29904
Human Swiss-Prot Number	O00418
Alternative Names	EEF2K; Eukaryotic elongation factor 2 kinase; eEF-2 kinase; eEF-2K; Calcium/calmodulin-dependent eukaryotic elongation factor 2 kinase