



ELK Biotechnology

GPR171 Rabbit pAb
Catalog Number: ES5698
For research use only.

Overview

Product Name	GPR171 polyclonal antibody
Host species	Rabbit
Applications	WB,IF,ELISA
Species Cross-Reactivity	Human,Mouse
Recommended dilutions	Western Blot:1/500-1/2000 Immunofluorescence:1/200-1/1000 ELISA:1/5000
Immunogen	Synthesized peptide derived from the Internal region of human GPR171. at AA range; 110-190
Specificity	GPR171 Polyclonal Antibody detects endogenous levels of GPR171 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	GPR171
Gene Name	GPR171
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	37kDa
Human Gene ID	29909
Human Swiss-Prot Number	O14626
Alternative Names	GPR171; H963; Probable G-protein coupled receptor 171; G-protein coupled receptor H963