



ELK Biotechnology

Histone H2B Rabbit pAb
Catalog Number: ES5715
For research use only.

Overview

Product Name	Histone H2B polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/40000
Immunogen	Synthesized peptide derived from human Histone H2B around the non-phosphorylation site of S14.
Specificity	Histone H2B Polyclonal Antibody detects endogenous levels of Histone H2B protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Histone H2B
Gene Name	HIST1H2BB
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	14kDa
Human Gene ID	3018
Human Swiss-Prot Number	P33778
Alternative Names	HIST1H2BB; H2BFF; Histone H2B type 1-B; Histone H2B.1; Histone H2B.f; H2B/f