

hnRNP C1/2 (phospho Ser260) rabbit pAb

Cat No.: ES5754

For research use only

Overview

Product Name hnRNP C1/2 (phospho Ser260) rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA **Species Cross-Reactivity** Human;Mouse;Rat

Recommended dilutions Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human hnRNP C1/2 around

the phosphorylation site of Ser260. AA

range:231-280

Specificity Phospho-hnRNP C1/2 (S260) Polyclonal Antibody

detects endogenous levels of hnRNP C1/2 protein

only when phosphorylated at S260.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Heterogeneous nuclear ribonucleoproteins C1/C2

Gene Name HNRNPC

Cellular localizationNucleus. Component of ribonucleosomes.PurificationThe antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 41kD
Human Gene ID 3183
Human Swiss-Prot Number P07910

Alternative Names HNRNPC; HNRPC; Heterogeneous nuclear

ribonucleoproteins C1/C2; hnRNP C1/C2

Background This gene belongs to the subfamily of ubiquitously

expressed heterogeneous nuclear


ribonucleoproteins (hnRNPs). The hnRNPs are RNA




+86-27-59760950 ELKbio@ELKbiotech.com wwv



binding proteins and they complex with heterogeneous nuclear RNA (hnRNA). These proteins are associated with pre-mRNAs in the nucleus and appear to influence pre-mRNA processing and other aspects of mRNA metabolism and transport. While all of the hnRNPs are present in the nucleus, some seem to shuttle between the nucleus and the cytoplasm. The hnRNP proteins have distinct nucleic acid binding properties. The protein encoded by this gene can act as a tetramer and is involved in the assembly of 40S hnRNP particles. Multiple transcript variants encoding at least two different isoforms have been described for this gene. [provided by RefSeq, Jul 2008],



Western Blot analysis of 293 cells using Phospho-hnRNP C1/2 (S260) Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003,Inventbiotech,MN,USA).




+86-27-59760950


Immunohistochemical analysis of paraffin-embedded Human brain. Antibody was diluted at 1:100(4° overnight). High-pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negetive contrl (right) obtaned from antibody was pre-absorbed by immunogen peptide.







Immunohistochemistry analysis of paraffin-embedded human brain, using hnRNP C1/2 (Phospho-Ser260) Antibody. The picture on the right is blocked with the phospho peptide.



+86-27-59760950

Western blot analysis of lysates from 293 cells treated with H2O2 100uM 15', using hnRNP C1/2 (Phospho-Ser260) Antibody. The lane on the right is blocked with the phospho peptide.

