



ELK Biotechnology

HSF1 Rabbit pAb
Catalog Number: ES5784
For research use only.

Overview

Product Name	HSF1 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/5000
Immunogen	Synthesized peptide derived from human HSF1 around the non-phosphorylation site of S307.
Specificity	HSF1 Polyclonal Antibody detects endogenous levels of HSF1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	HSF1
Gene Name	HSF1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	30kDa
Human Gene ID	3297
Human Swiss-Prot Number	Q00613
Alternative Names	HSF1; HSTF1; Heat shock factor protein 1; HSF 1; Heat shock transcription factor 1; HSTF 1