



ELK Biotechnology

Lambda 5 Rabbit pAb
Catalog Number: ES5863
For research use only.

Overview

Product Name	Lambda 5 polyclonal antibody
Host species	Rabbit
Applications	IHC-p,WB,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/40000
Immunogen	Synthesized peptide derived from the C-terminal region of human Lambda 5. at AA range; 130-210
Specificity	Lambda 5 Polyclonal Antibody detects endogenous levels of Lambda 5 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Lambda 5
Gene Name	IGLL1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	3543
Human Swiss-Prot Number	P15814
Alternative Names	IGLL1; IGL1; Immunoglobulin lambda-like polypeptide 1; CD179 antigen-like family member B; Ig lambda-5; Immunoglobulin omega polypeptide; Immunoglobulin-related protein 14.1; CD antigen CD179b