

MRLC2 (phospho Ser20) rabbit pAb

Cat No.:ES5880

For research use only

Overview

Product Name	MRLC2 (phospho Ser20) rabbit pAb
Host species	Rabbit
Applications	WB;IHC;IF;ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	Western Blot: 1/500 - 1/2000.
	Immunohistochemistry: 1/100 - 1/300.
	Immunofluorescence: 1/200 - 1/1000. ELISA:
	1/20000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized
	peptide derived from human Myosin regulatory light
	chain 2 around the phosphorylation site of Ser18.
	AA range:3-52
Specificity	Phospho-MRLC2 (S20) Polyclonal Antibody detects
	endogenous levels of MRLC2 protein only when
	phosphorylated at S20.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Myosin regulatory light polypeptide 9
Gene Name	MYL9
Cellular localization	Cytoplasm, cytoskeleton . Cytoplasm, cell cortex .
	Colocalizes with F-actin, MYH9 and PIEZO1 at the
	actomyosin cortex in myoblasts.
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
Closelity	epitope-specific immunogen.
Clonality Concentration	Polyclonal
Observed band	1 mg/ml 18kD
Human Gene ID	10398/10627
Human Swiss-Prot Number	P24844
Alternative Names	MYL9; MLC2; MRLC1; MYRL2; Myosin regulatory
AILEI HALIVE INAILLES	light polypeptide 9; 20 kDa myosin light chain; LC20;
	ingine poryperpende 5, 20 kba myosin ingine enam, 1020,



+86-27-59760950

ELKbio@ELKbiotech.com


www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C




Background

MLC-2C; Myosin RLC; Myosin regulatory light chain 2; smooth muscle isoform; Myosin regulatory light chain 9; Myosin regulatory light chain MRL Myosin, a structural component of muscle, consists of two heavy chains and four light chains. The protein encoded by this gene is a myosin light chain that may regulate muscle contraction by modulating the ATPase activity of myosin heads. The encoded protein binds calcium and is activated by myosin light chain kinase. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008],



Enzyme-Linked Immunosorbent Assay (Phospho-ELISA) for Immunogen Phosphopeptide (Phospho-left) and Non-Phosphopeptide (Phospho-right), using Myosin regulatory light chain 2 (Phospho-Ser18) Antibody



Immunofluorescence analysis of HUVEC cells, using Myosin regulatory light chain 2 (Phospho-Ser18) Antibody. The picture on the right is blocked with the phospho peptide.



+86-27-59760950


ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



Myosin regulatory light chain 2 (pSer18)-- Immunohistochemistry analysis of paraffin-embedded human brain, using Myosin regulatory light chain 2 (Phospho-Ser18) Antibody. The picture on the right is blocked with the phospho peptide.



blocked with the phospho peptide.



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C