

IL-8Rβ Rabbit pAb Catalog Number: ES5887 For research use only.

Overview

Product Name	IL-8Rβ polyclonal antibody
Host species	Rabbit
Applications	IF,WB,ELISA
Species Cross-Reactivity	Human,Mouse
Recommended dilutions	Western Blot:1/500-1/2000 Immunofluorescence:1/200-1/1000 ELISA:1/10000
Immunogen	Synthesized peptide derived from human IL-8R β around the non-phosphorylation site of S347.
Specificity	IL-8R β Polyclonal Antibody detects endogenous levels of IL-8R β protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	IL-8Rβ
Gene Name	CXCR2
Purification	The antibody was affinity-purified from rabbit antiserum by affinity- chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	40kDa
Human Gene ID	3579
Human Swiss-Prot Number	P25025
Alternative Names	CXCR2; IL8RB; C-X-C chemokine receptor type 2; CXC-R2; CXCR-2; CDw128b; GRO/MGSA receptor; High affinity interleukin-8 receptor B; IL-8R B; IL-8 receptor type 2; CD antigen CD182