



# ELK Biotechnology

**IL-9R Rabbit pAb**  
**Catalog Number: ES5890**  
**For research use only.**

## Overview

---

Product Name	IL-9R polyclonal antibody
Host species	Rabbit
Applications	IHC-p,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	<b>Immunohistochemistry:1/100-1/300</b> <b>ELISA:1/20000</b>
Immunogen	Synthesized peptide derived from human IL-9R around the non-phosphorylation site of S519.
Specificity	IL-9R Polyclonal Antibody detects endogenous levels of IL-9R protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	IL-9R
Gene Name	IL9R
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	<b>kDa</b>
Human Gene ID	3581
Human Swiss-Prot Number	Q01113
Alternative Names	IL9R; Interleukin-9 receptor; IL-9 receptor; IL-9R; CD antigen CD129