



## ELK Biotechnology

### IRF-3 Rabbit pAb

Catalog Number: ES5925

For research use only.

#### Overview

---

Product Name	IRF-3 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/20000
Immunogen	Synthesized peptide derived from human IRF-3 around the non-phosphorylation site of S386.
Specificity	IRF-3 Polyclonal Antibody detects endogenous levels of IRF-3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	IRF-3
Gene Name	IRF3
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	50kDa
Human Gene ID	3661
Human Swiss-Prot Number	Q14653
Alternative Names	IRF3; Interferon regulatory factor 3; IRF-3