



ELK Biotechnology

IRS-1 Rabbit pAb
Catalog Number: ES5932
For research use only.

Overview

Product Name	IRS-1 polyclonal antibody
Host species	Rabbit
Applications	IHC-p,IF,WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 Immunofluorescence:1/200-1/1000 ELISA:1/20000
Immunogen	Synthesized peptide derived from human IRS-1 around the non-phosphorylation site of S307.
Specificity	IRS-1 Polyclonal Antibody detects endogenous levels of IRS-1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	IRS-1
Gene Name	IRS1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	3667
Human Swiss-Prot Number	P35568
Alternative Names	IRS1; Insulin receptor substrate 1; IRS-1