



ELK Biotechnology

AR Rabbit pAb
Catalog Number: ES5939
For research use only.

Overview

Product Name	AR polyclonal antibody
Host species	Rabbit
Applications	IHC-p,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Immunohistochemistry:1/100-1/300 ELISA:1/40000
Immunogen	Synthesized peptide derived from human AR around the non-phosphorylation site of S213.
Specificity	AR Polyclonal Antibody detects endogenous levels of AR protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	AR
Gene Name	AR
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	367
Human Swiss-Prot Number	P10275
Alternative Names	AR; DHTR; NR3C4; Androgen receptor; Dihydrotestosterone receptor; Nuclear receptor subfamily 3 group C member 4