



ELK Biotechnology

Cytokeratin 14 Rabbit pAb
Catalog Number: ES6032
For research use only.

Overview

Product Name	Cytokeratin 14 polyclonal antibody
Host species	Rabbit
Applications	IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Immunohistochemistry:1/100-1/300 ELISA:1/10000
Immunogen	Synthesized peptide derived from the Internal region of human Cytokeratin 14.
Specificity	Cytokeratin 14 Polyclonal Antibody detects endogenous levels of Cytokeratin 14 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Cytokeratin 14
Gene Name	KRT14
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	3861
Human Swiss-Prot Number	P02533
Alternative Names	KRT14; Keratin; type I cytoskeletal 14; Cytokeratin-14; CK-14; Keratin-14; K14