



ELK Biotechnology

GPR153 Rabbit pAb
Catalog Number: ES6037
For research use only.

Overview

Product Name	GPR153 polyclonal antibody
Host species	Rabbit
Applications	WB,IF,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 Immunofluorescence:1/200-1/1000 ELISA:1/5000
Immunogen	Synthesized peptide derived from the Internal region of human GPR153. at AA range; 260-340
Specificity	GPR153 Polyclonal Antibody detects endogenous levels of GPR153 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	GPR153
Gene Name	GPR153
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	65kDa
Human Gene ID	387509
Human Swiss-Prot Number	Q6NV75
Alternative Names	GPR153; PGR1; Probable G-protein coupled receptor 153; G-protein coupled receptor PGR1