

Rho G rabbit pAb

Cat No.: ES6076

For research use only

Overview

Product Name Rho G rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA Species Cross-Reactivity Human;Mouse;Rat

Recommended dilutions Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human RHOG. AA

range:97-146

Specificity Rho G Polyclonal Antibody detects endogenous

levels of Rho G protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Rho-related GTP-binding protein RhoG

Gene Name RHOG

Cellular localizationCell membrane ; Lipid-anchor ; Cytoplasmic side .PurificationThe antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 23kD
Human Gene ID 391
Human Swiss-Prot Number P84095

Alternative Names RHOG; ARHG; Rho-related GTP-binding protein

RhoG

Background This gene encodes a member of the Rho family of

small GTPases, which cycle between inactive GDP-bound and active GTP-bound states and function as molecular switches in signal


transduction cascades. Rho proteins promote

+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com



reorganization of the actin cytoskeleton and regulate cell shape, attachment, and motility. The encoded protein facilitates translocation of a functional guanine nucleotide exchange factor (GEF) complex from the cytoplasm to the plasma membrane where ras-related C3 botulinum toxin substrate 1 is activated to promote lamellipodium formation and cell migration. Two related pseudogene have been identified on chromosomes 20 and X. [provided by RefSeq, Aug 2011],



Western blot analysis of lysates from Jurkat cells, using RHOG Antibody. The lane on the right is blocked with the synthesized peptide.



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).

