



ELK Biotechnology

CD107a Rabbit pAb
Catalog Number: ES6085
For research use only.

Overview

Product Name	CD107a polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/40000
Immunogen	Synthesized peptide derived from the Internal region of human CD107a. at AA range; 90-170
Specificity	CD107a Polyclonal Antibody detects endogenous levels of CD107a protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	CD107a
Gene Name	LAMP1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	125kDa
Human Gene ID	3916
Human Swiss-Prot Number	P11279
Alternative Names	LAMP1; Lysosome-associated membrane glycoprotein 1; LAMP-1; Lysosome-associated membrane protein 1; CD107 antigen-like family member A; CD antigen CD107a