



ELK Biotechnology

Myoglobin Rabbit pAb Catalog Number: ES6198 For research use only.

Overview

Product Name	Myoglobin polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/5000
Immunogen	Synthesized peptide derived from the Internal region of human Myoglobin. at AA range; 40-120
Specificity	Myoglobin Polyclonal Antibody detects endogenous levels of Myoglobin protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Myoglobin
Gene Name	MB
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	17kDa
Human Gene ID	4151
Human Swiss-Prot Number	P02144
Alternative Names	MB; Myoglobin