



ELK Biotechnology

Mcl-1 Rabbit pAb
Catalog Number: ES6206
For research use only.

Overview

Product Name	Mcl-1 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Monkey
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/10000
Immunogen	Synthesized peptide derived from human Mcl-1 around the non-phosphorylation site of S159.
Specificity	Mcl-1 Polyclonal Antibody detects endogenous levels of Mcl-1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Mcl-1
Gene Name	MCL1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	35kDa
Human Gene ID	4170
Human Swiss-Prot Number	Q07820
Alternative Names	MCL1; BCL2L3; Induced myeloid leukemia cell differentiation protein Mcl-1; Bcl-2-like protein 3; Bcl2-L-3; Bcl-2-related protein EAT/mcl1; mcl1/EAT