



ELK Biotechnology

MEK Kinase-1 Rabbit pAb
Catalog Number: ES6229
For research use only.

Overview

Product Name	MEK Kinase-1 polyclonal antibody
Host species	Rabbit
Applications	WB,IF,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 Immunofluorescence:1/200-1/1000 ELISA:1/5000
Immunogen	Synthesized peptide derived from the N-terminal region of human MEK Kinase-1. at AA range; 150-230
Specificity	MEK Kinase-1 Polyclonal Antibody detects endogenous levels of MEK Kinase-1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	MEK Kinase-1
Gene Name	MAP3K1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	165kDa
Human Gene ID	4214
Human Swiss-Prot Number	Q13233
Alternative Names	MAP3K1; MAPKKK1; MEKK; MEKK1; Mitogen-activated protein kinase kinase kinase 1; MAPK/ERK kinase kinase 1; MEK kinase 1; MEKK 1