



ELK Biotechnology

c-Myb Rabbit pAb
Catalog Number: ES6299
For research use only.

Overview

Product Name	c-Myb polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Cross-Reactivity	Human,Mouse
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 Immunofluorescence:1/200-1/1000 ELISA:1/10000
Immunogen	Synthesized peptide derived from human c-Myb around the non-phosphorylation site of S12.
Specificity	c-Myb Polyclonal Antibody detects endogenous levels of c-Myb protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	c-Myb
Gene Name	MYB
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	72kDa
Human Gene ID	4602
Human Swiss-Prot Number	P10242
Alternative Names	MYB; Transcriptional activator Myb; Proto-oncogene c-Myb