



ELK Biotechnology

MLTK Rabbit pAb
Catalog Number: ES6566
For research use only.

Overview

Product Name	MLTK polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/40000
Immunogen	Synthesized peptide derived from human MLTK around the non-phosphorylation site of S165.
Specificity	MLTK Polyclonal Antibody detects endogenous levels of MLTK protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	MLTK
Gene Name	MLTK
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	91kDa
Human Gene ID	51776
Human Swiss-Prot Number	Q9NYL2
Alternative Names	MLTK; ZAK; HCCS4; Mitogen-activated protein kinase kinase kinase MLT; Human cervical cancer suppressor gene 4 protein; HCCS-4; Leucine zipper-and sterile alpha motif-containing kinase; MLK-like mitogen-activated protein triple kinase; Mixe