



ELK Biotechnology

PR Rabbit pAb
Catalog Number: ES6582
For research use only.

Overview

Product Name	PR polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Monkey
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/5000
Immunogen	Synthesized peptide derived from human PR around the non-phosphorylation site of S190.
Specificity	PR Polyclonal Antibody detects endogenous levels of PR protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	PR
Gene Name	PGR
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	99kDa
Human Gene ID	5241
Human Swiss-Prot Number	P06401
Alternative Names	PGR; NR3C3; Progesterone receptor; PR; Nuclear receptor subfamily 3 group C member 3