



# ELK Biotechnology

**PI 3-kinase p85 $\beta$  Rabbit pAb**  
**Catalog Number: ES6593**  
**For research use only.**

## Overview

---

<b>Product Name</b>	PI 3-kinase p85 $\beta$ polyclonal antibody
<b>Host species</b>	Rabbit
<b>Applications</b>	<b>IHC-p,ELISA</b>
<b>Species Cross-Reactivity</b>	<b>Human,Mouse,Rat</b>
<b>Recommended dilutions</b>	<b>Immunohistochemistry:1/100-1/300</b> <b>ELISA:1/10000</b>
<b>Immunogen</b>	Synthesized peptide derived from the Internal region of human PI 3-kinase p85 $\beta$ .
<b>Specificity</b>	PI 3-kinase p85 $\beta$ Polyclonal Antibody detects endogenous levels of PI 3-kinase p85 $\beta$ protein.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Protein Name</b>	PI 3-kinase p85 $\beta$
<b>Gene Name</b>	PIK3R2
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	<b>kDa</b>
<b>Human Gene ID</b>	5296
<b>Human Swiss-Prot Number</b>	O00459
<b>Alternative Names</b>	PIK3R2; Phosphatidylinositol 3-kinase regulatory subunit beta; PI3-kinase regulatory subunit beta; PI3K regulatory subunit beta; PtdIns-3-kinase regulatory subunit beta; Phosphatidylinositol 3-kinase 85 kDa regulatory subunit beta; PI3-kinase