



ELK Biotechnology

PKD2 Rabbit pAb
Catalog Number: ES6598
For research use only.

Overview

Product Name	PKD2 polyclonal antibody
Host species	Rabbit
Applications	IF,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Immunofluorescence:1/200-1/1000 ELISA:1/10000
Immunogen	Synthesized peptide derived from human PKD2 around the non-phosphorylation site of S812.
Specificity	PKD2 Polyclonal Antibody detects endogenous levels of PKD2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	PKD2
Gene Name	PKD2
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	5311
Human Swiss-Prot Number	Q13563
Alternative Names	PKD2; Polycystin-2; Autosomal dominant polycystic kidney disease type II protein; Polycystic kidney disease 2 protein; Polycystin; R48321