

EDG-8 rabbit pAb

Cat No.:ES6624

For research use only

Overview

Product Name EDG-8 rabbit pAb

Host species Rabbit
Applications IHC;IF;ELISA

Species Cross-Reactivity Human; Rat; Mouse;

Recommended dilutions Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human EDG8. AA

range:335-384

Specificity EDG-8 Polyclonal Antibody detects endogenous

levels of EDG-8 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Sphingosine 1-phosphate receptor 5

Gene Name S1PR5

Cellular localizationCell membrane; Multi-pass membrane protein.PurificationThe antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

Observed band

Human Gene ID 53637 Human Swiss-Prot Number Q9H228

Alternative Names S1PR5; EDG8; Sphingosine 1-phosphate receptor 5;

S1P receptor 5; S1P5; Endothelial differentiation G-protein-coupled receptor 8; Sphingosine

1-phosphate receptor Edg-8; S1P receptor Edg-8

Background The lysosphingolipid sphingosine 1-phosphate (S1P)

regulates cell proliferation, apoptosis, motility, and

neurite retraction. Its actions may be both




+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.c




intracellular as a second messenger and extracellular as a receptor ligand. S1P and the structurally related lysolipid mediator lysophosphatidic acid (LPA) signal cells through a set of G protein-coupled receptors known as EDG receptors. Some EDG receptors (e.g., EDG1; MIM 601974) are S1P receptors; others (e.g., EDG2; MIM 602282) are LPA receptors. [supplied by OMIM, Mar 2008],

Immunofluorescence analysis of A549 cells, using EDG8 Antibody. The picture on the right is blocked with the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using EDG8 Antibody. The picture on the right is blocked with the synthesized peptide.



+86-27-59760950

