



ELK Biotechnology

Histone H2B Rabbit pAb
Catalog Number: ES6641
For research use only.

Overview

Product Name	Histone H2B polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/10000
Immunogen	Synthesized peptide derived from human Histone H2B around the non-acetylation site of K15. at AA range; 1-80
Specificity	Histone H2B Polyclonal Antibody detects endogenous levels of Histone H2B protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Histone H2B
Gene Name	H2BFS
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	20kDa
Human Gene ID	54145
Human Swiss-Prot Number	P57053
Alternative Names	H2BFS; Histone H2B type F-S; Histone H2B.s; H2B/s