

NLE1 rabbit pAb

Cat No.: ES6655

For research use only

Overview

Product Name NLE1 rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA **Species Cross-Reactivity** Human;Monkey

Recommended dilutions Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human NLE1. AA

range:301-350

Specificity NLE1 Polyclonal Antibody detects endogenous levels

of NLE1 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Notchless protein homolog 1

Gene Name NLE1

Cellular localization Nucleus, nucleolus.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

ClonalityPolyclonalConcentration1 mg/mlObserved band48kDHuman Gene ID54475Human Swiss-Prot NumberQ9NVX2

+86-27-59760950

Alternative Names NLE1; HUSSY-07; Notchless protein homolog 1


Background similarity:Contains 8 WD repeats.,




ELKbio@ELKbiotech.com

www.elkbiotech.com





Western blot analysis of lysates from COS cells, using NLE1 Antibody. The lane on the right is blocked with the synthesized peptide.



Immunohistochemical analysis of paraffin-embedded human brain. 1, Tris-EDTA,pH9.0 was used for antigen retrieval. 2 Antibody was diluted at 1:200(4° overnight.3,Secondary antibody was diluted at 1:200(room temperature, 45min).

